



VK Healthcare System

**Consumable And
Accessories**

Brochure

Insert Earphones IP 30

----- X

The IP30 Insert phone from RadioEar has been designed to meet the ISO 389-2 standard and matches the performance of the most popular insert phones used in the industry today. These devices have proved to be quite accurate in a variety of testing conditions and are very effective in filtering out unwanted background noise. This makes them an ideal choice for nonideal examination settings. Additionally the RadioEar IP30 Insert Phone does not require the use of bulky uncomfortable headbands.



DD 45 Headphones

----- X

- Engineered to match the performance and function of the most popular earphones available for audiological testing in today's marketplace.
- The DD45 air transducer is manufactured to the highest quality standards and is a more economical alternative.
- RadioEar has an extensive history of successfully serving audiometer manufacturers and service providers. We recognize the critical importance these devices have for our customers.
- The DD45 Audiometric Earphones and all of the accessories are manufactured for RadioEar.
- These earphones are fully supported by RadioEar, and receive the same focused attention and personal service that we provide our customers for bone transducers, headbands and cords.



B71 Bone Conduction

----- X

Since its introduction in the 1970s the RadioEar B71 has been the bone transducer of choice for audiometric testing. As a consequence of ever changing standards we continuously endeavor to update our product, so that it is compliant with all the latest regulations. RadioEar has an extensive history of successfully serving the audiometric industry. We recognize the critical importance these devices have for our customers, and are sincerely dedicated to serving them.



Nuprep Gel

----- X



Weaver Nuprep Skin Preparation Gel is easy and convenient to use. A flip-top cap lets you dispense the gel with one hand, leaving the other free. The flat top lets you stand the tube on its head within easy reach. Nuprep gel has a consistent texture, so it won't

blurt when applied. You can dispense precisely the amount you need. Nuprep does not contain Latex in either the product ingredients or packaging, and was manufactured using latex-free gloves.

Ten-20 Paste

----- X

Weaver Ten20 Conductive EEG Paste will hold the electrodes in place so there is no need to disrupt a test to reattach the electrodes. This reliable paste will conduct the signal that the technologist/technician is relying on for clear results.

Expert technologists and technicians throughout the world rely on the optimal





balance of adhesiveness and conductivity of Ten20 paste. Apply the electrode using Ten20, and it will remain in place until it is removed.

Disposable Sanibel™ ADI ear tips Box

-----X

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip:

Amplivox: all Otowave versions.

Interacoustics: Titan, AT235, AA220, Eclipse, MT10, Sera, OtoRead II. **MAICO:** easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, EroScan Plus. Alternative to AKR (KR) / IA series: Old MAICO MI-series + RaceCar black probe versions together with adaptor kit (8500430). **GSI:** Novus, TympStar Pro, Corti and Allegro. **Natus:** Zodiac (not 901) and AccuScreen.



Probe Tips

-----X

Probe tip for OtoWave™ and MT10™	
Probe tip for Titan™ and easyTymp™	



Probe tip Kit	
ER10D probe tips for ERO•SCAN™ Pro and OtoRead™	
Probe tip kit for AT235 Old and MI 34	
Probe Tube for MAICO Ero Scan Plus, OtoRead, GSI Corti	

Foam Ear Tip 10mm and 13 mm

----- x

Sanibel insert foam tips expand quickly to ensure a great fit. Once inserted into the ear canal, the expanded foam ensures a sturdy fit that keeps the insert phone in place during the entire measurement.

- Used for AUD/IMP/ABR
- Compatible with all types of insert phones
- Conforms to Medical Device Regulations
- Latex Free



- Single patient use only

Transparent tube / Sound tubes for insert earphone

----- X

Sound tubes are compatible with all ABR equipment. Among the many benefits these sound tubes provide are improved ambient noise attenuation, facilitated infection control, and improved audiometry test and retest reliability.



- Replacement sound tubes with nipple connectors
- 9.85" front tube
- Polyurethane

Gold Eartip Electrode 10mm and 13mm for ECoG

----- X

Gold foil tiptrodes serve as electrodes to deliver stimulus via the silicone tubing connected to the insert earphone. Optimally recorded evoked potentials from the ear canal and are disposable after use.



- 13mm and 10 mm
- Disposable electrode combined with an insert earphone foam tip
- Used to record evoked potentials from the ear canal
- Standard size
- TTE25 Gold foil wrapped foam

Button electrode cable

-----X

Imported Button Electrode Cable available in different colors with good impedance matching, used for BERA, ECG, EKG, Holter.



Jumper cable

-----X


Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.



Gold Cup Ep Electrode (1.5M)

-----X

Gold cup electrodes on a color-coded, pull-apart ribbon cable. Use conductive paste to collect EEG, EMG, or ECG data. Compatible with ALL OpenBCI boards.



Probe Tip cleaner

-----X

This Bridge & Implant floss is used for cleaning the Probe tip kit of OAE, Bera and Impedance Audiometer.





Impression Syringe

-----X

Easy to fill, inject and clean. Provides a stable and easy injection. 4 mm. Our friends at Pacific Audiology Group have launched a new in-depth online course; "Comprehensive Earmold Impressions".

The contents are relevant to every Hearing Healthcare Professional who wishes to receive competent training and expert knowledge of the entire process, from start to finish, with all aspects of impression-taking.



Replacement cable for old version of B-71

-----X



DD45™/TDH39™ replacement cable, symmetric

-----X

Compatible with: Benson CCA-100mini, Benson CCA-200mini, GSI 61, GSI AudioStar, GSI AudioStar Pro, GSI GSI38, GSI GSI39, GSI Pello, Interacoustics AC40, Interacoustics AD226 New, Interacoustics AD528, Interacoustics AD629, Interacoustics Affinity Compact, Interacoustics Callisto, Interacoustics Equinox, MAICO MI36, Natus Madsen Aurical, Smart Diagnostic Devices Smart Tone, MedRx AVANT A2D+.





Thermal Paper

- - - - - X

Paper for thermo printer. Fits GSI33/38, Madsen Zodiac 901, Interacoustics AT235 / AA220 / AA222 / AZ26, ATA22t / AC40 and Maico easyTymp, MI22, MI23, MI24, MI26, MI34 and MI 44 and many more. Fits most popular audiometers and tympanometers.



For further queries Please contact us :

V K Healthcare System

Contact Number: 7678531115

Email: vkhealthcaresystem@gmail.com

Website: www.vkhealthcaresystem.com

